



QUART dent/digitest

Dental X-Ray QA/QC Test Phantom



QUART dent/digitest

Test Phantom for Quality Assurance and Critical Examination of Dental X-Ray Equipment



Only one exposure is required to collect all necessary parameters to determine the imaging quality of the x-ray equipment. After the exposure, the test image is evaluated visually.

The *QUART dent/digitest* line of phantoms fully comply with IEC 61223-3-4 and IEC 61223-2-7, DIN 6868-151 and DIN 6868-5 as well as IPEM, ÖNORM, PN-EN standards and AFSSAPS regulations.

Phantom and Test Procedure

QC tests are carried out at initial equipment installation to establish a visual reference. On a regular routine basis, they are meant to ensure the constancy of the system's imaging capabilities. Any deterioration in imaging performance would be revealed by the visual evaluation of the test images. A single image of the QA/QC phantom contains a lot of crucial information for that purpose.

Test Parameters

- Spatial Resolution
- _ Low-Contrast Resolution
- _ Radiation Field Alignment
- _ Image Homogeneity
- _ Dose Equivalency Check
- _ Artefacts, Image Flaws, etc.

Customised Holder Design

For QA/QC performance at panoramic x-ray equipment a customised holder can be provided for fast and reproducible positioning of the phantom in the tomographic plane. Adaptors for a variety of systems are available:

AJAT ART Plus Series
BlueX PantOs Series

GENDEX Orthoralix 8500/9200

J. MORITA IC-5/Veraviewepocs

KODAK-CARESTREAM K8000/K9000

MYRAY Hyperion X Series
OWANDY Imax Series

PLANMECA Proline/ProOne/ProMax

RAY MEDICAL Rayscan alpha Series


SIRONA OP100/OP200/Orthophos Series

SOREDEX OP30/OP100/OP200
VATECH/EWOO Pax/Picasso Series
VILLA Sistemi Rotograph Series

NB: The holder can be customised for any panoramic system.









QUART dent/digitest

Test Phantom for Quality Assurance and Critical Examination of Dental X-Ray Equipment

A Genuine History of Invention

In the mid 1990's QUART developed and designed the first dent/digitest QA/QC phantom. This phantom soon became the National Standard in Germany when it was described and represented in the DIN standard text. Some time later it was also discussed by the European IEC standard committee and became a standardised solution for Europe as well.

Until today the phantom has not changed much in its appearance. However, QUART added genuine features to further enhance the ease of use of the phantom: the separable design for precise panoramic and cephalometric testing, an enhanced geometric test, a customisable panoramic holder and a special holder solution for cephalometric systems.

Phantom Description

- _ Separable Design for all dental modalities
- _ Size: 80x80x36 mm (LxWxH)
- _ Slit to slide out intra-oral sensor without pulling the cable; also for secure sensor cable routing
- _ 2-point fitting for cephalometric systems
- _ V-nut for optional dent/digi H OPG holder fit (see accessories below)
- _ 3 centring rings with different diameters to fit standard tubes
- _ Slots for digital storage screen, intra-oral sensor, dose detector
- _ 45° Line Pair test (2.5/2.8/3.1/5.0/5.8/6.3 Lp/mm; 0.05 Pb)
- _ 4 Low Contrast objects (Ø 2.5/2/1.5/1 mm)
- _ Integrated filtration 6 mm AI (purity 99.5 % guaranteed)

Delivery includes

- ✓ QUART dent/digitest 2.1 or 3.1 Test Phantom
- ✓ QUART dent/digi H Phantom Holder (optional)
- ✓ QUART cu08 Added Filtration for pano/ceph (optional)
- ✓ Manuals
- ✓ Transport Case with Foam Insert

Also available

Variety of customised Holders

QUART dent/digi H uni

Universal tripod stand

QUART dentII Phantom for Conventional Dental QA/QC

QUART dentII mpc Phantom Holder for conventional pano/ceph

QUART cu10d Added DIN complient filtration
QUART dent/digitest M1/M2 Manufacturer Test Phantom
QUART dido/EASY Constancy Test Dosemeter









www.quart.de

